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XIEM
EX

Fusion

When a new idea meets
ZEOCLIP® technology

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
Figure1


Figure2

1. Sakamoto N, Osada T, Shibuya T, et al. The facilitation of a new traction device (S-O clip) assisting endoscopic submucosal dissection for superficial colorectal neoplasms. Endoscopy. 2008 Sep;40 Suppl 2:E94-5.
 2. Sakamoto N, Osada T, Shibuya T, et al. Endoscopic submucosal dissection of large colorectal tumors by using a novel spring-action S-O clip for traction (with video). Gastrointest Endosc. 2009 Jun;69(7):1370-4.
 3.Ritsuno H, Sakamoto N, Osada T,et al. Prospective clinical trial of traction device-assisted endoscopic submucosal dissection of large superficial colorectal tumors using the S-O clip.Surgical Endoscopy 2014 in press.
 4.Takeda T, Murakami T, Sakamoto N,et al. Traction device to remove an adenoma in the appendiceal orifice by endoscopic submucosal dissection. Endoscopy. 2013; 45 :E239-40.
 5.Sakamoto N, Beppu K, Matsumoto K, et al. "Loop Clip", a new closure device for large mucosal defects after EMR and ESD. Endoscopy. 2008;40:E97-8.
 6.Mori H, Sakamoto N, Osada T, et al. The 'Loop Clip' is useful for closing large mucosal defects after colorectal endoscopic submucosal dissection: a preliminary clinical study. Dig Endosc. 2011 Oct; 23(4):330-1.

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1. Ono H, Kondo H, Gotoda T et al. Endoscopic mucosal resection for treatment of early gastric cancer. *Gut* 2001 ; 48 : 225-9.
2. Yamamoto H, Kawata H, Sunada K et al. Successful en-bloc resection of large superficial tumors in the stomach and colon using sodium hyaluronate and small-caliber-tip transparent hood. *Endoscopy* 2003 ; 35 : 690-4.
3. Oyama T, Tomori A, Hotta K et al. Endoscopic submucosal dissection of early esophageal cancer. *Clin Gastroenterol Hepatol* 2005 ; 3 (7 Suppl 1) : S67-S70.
4. Saito Y, Emura F, Matsuda T et al. A new sinker-assisted endoscopic submucosal dissection for colorectal tumors. *Gastrointest Endosc* 2005 ; 62 : 297-301.
5. Yamamoto K, Hayashi S, Saiki H et al. Endoscopic submucosal dissection for large superficial colorectal tumors using the "clip-flap method". *Endoscopy* 2015 ; 47 : 262-5.
6. Imaeda H, Iwao Y, Ogata H et al. A new technique for endoscopic submucosal dissection for early gastric cancer using an external grasping forceps. *Endoscopy* 2006 ; 38 : 1007-10.
7. Teoh AY, Chiu PW, Hon SF et al. Ex vivo comparative study using the Endolifter® as a traction device for enhancing submucosal visualization during endoscopic submucosal dissection. *Surg Endosc* 2013 ; 27 : 1422-7.
8. Oyama T. Counter traction makes endoscopic submucosal dissection easier. *Clin Endosc* 2012 ; 45 : 375-8.
9. Gotoda T, Oda I, Tamakawa K et al. Prospective clinical trial of magnetic-anchor-guided endoscopic submucosal dissection for large early gastric cancer (with videos). *Gastrointest Endosc* 2009 ; 69 : 10-5.
10. Yonezawa J, Kaise M, Sumiyama K et al. A novel double-channel therapeutic endoscope ("R-scope") facilitates endoscopic submucosal dissection of superficial gastric neoplasms. *Endoscopy* 2006 ; 38 : 1011-5.
11. Neuhaus H, Costamagna G, Devière J et al. Endoscopic submucosal dissection (ESD) of early neoplastic gastric lesions using a new double-channel endoscope (the "R-scope"). *Endoscopy* 2006 ; 38 : 1016-23.
12. Uraoka T, Kato J, Ishikawa S et al. Thin endoscope-assisted endoscopic submucosal dissection for large colorectal tumors (with videos). *Gastrointest Endosc* 2007 ; 66 : 836-9.
13. Sakamoto N, Osada T, Shibuya T et al. The facilitation of a new traction device (S-O clip) assisting endoscopic submucosal dissection for superficial colorectal neoplasms. *Endoscopy* 2008 ; 40 : E94-5.
14. Sakamoto N, Osada T, Shibuya T et al. Endoscopic submucosal dissection of large colorectal tumors by using a novel spring-action S-O clip for traction (with video). *Gastrointest Endosc* 2009 ; 69 : 1370-4.
15. Ota M, Nakamura T, Hayashi K et al. Usefulness of clip traction in the early phase of esophageal endoscopic submucosal dissection. *Dig Endosc* 2012 ; 24 : 315-8.
16. Takeda T, Murakami T, Sakamoto N et al. Traction device to remove an adenoma in the appendiceal orifice by endoscopic submucosal dissection. *Endoscopy* 2013 ; 45 : E239-40.
17. Ritsuno H, Sakamoto N, Osada T et al. Large superficial tumor of the colon involving a diverticulum removed by endoscopic submucosal dissection. *Gastrointest Endosc* 2015 ; 82 : 751.
18. Ritsuno H, Sakamoto N, Osada T et al. Prospective clinical trial of traction device-assisted endoscopic submucosal dissection of large superficial colorectal tumors using the S-O clip. *Surg Endosc* 2014 ; 28 : 3143-9.

ENDOSCOPIC SUBMUCOSAL DISSECTION OF COLORECTAL TUMORS USING A TRACTION DEVICE CALLED THE S-O CLIP

Naoto SAKAMOTO¹⁾, Taro OSADA³⁾, Hideaki RITSUNO¹⁾,
Takashi MURAKAMI¹⁾, Tomoyoshi SHIBUYA¹⁾, Yuichi TOMIKI²⁾
AND Sumio WATANABE¹⁾

1) *Department of Gastroenterology, Juntendo University School of Medicine.*

2) *Department of Coloproctological Surgery, Juntendo University School of Medicine.*

3) *Department of Gastroenterology, Juntendo University Urayasu Hospital.*

Endoscopic submucosal dissection (ESD) allows en bloc resection of large gastrointestinal neoplasms regardless of their size. However, ESD is difficult to perform and it is associated with a longer procedure time and higher risk of complications such as perforation than endoscopic mucosal resection (EMR). To address this difficulty, several traction methods have been introduced to facilitate ESD procedure.

We designed a traction device called the S-O clip (Zeon Medical Co., Ltd., Tokyo, Japan) and reported it previously. The S-O clip consists of a metallic clip attached to the end of a spring, which is connected at its other end to a nylon loop. It can be passed through an endoscope instrument channel and can easily pull a lesion at any location up without withdrawing the endoscope. We describe the basic and recommended usage of the S-O clip.

1. A circumferential incision around the tumor in the submucosal layer is performed.
2. After separation of the tumor from the surrounding normal mucosa, an S-O clip is attached to the edge of the exfoliated mucosa.
3. A regular clip is used to grasp the distal nylon loop attached to the S-O clip and then the regular clip is inserted and applied to the colon wall opposite the lesion, enabling traction and therefore opening the resection margin.
4. The traction applied by the S-O clip pulls on the edge of the lesion, allowing the cutting line in the submucosal layer to be visualized and resulting in safe and successful en bloc dissection.
5. After dissection, the nylon loop is cut with a knife, and the specimen is brought out using the forceps. The endoscope is not withdrawn during the entire procedure.

S-O clip-assisted ESD is safe and effective for en bloc resection of large superficial neoplasms in the colorectum.